CORTEX TECHNOLOGY DENMARK

DE: Sprühhebel
 DK: Triggerhåndtag
 EE: Käepide
 ES: Pulsador
 FI: Käsi kahva
 FR: Gâchette
 GR: Ααβή σκανόαλης
 HU: Indító kar
 IT: Leva
 NL: Trekhandvat
 NO: Aktiveringshandtak
 PL: Spust ręczny
 PT: Disparador
 SE: Ventil handtag
 RO: Mâner de declanșare
 UK: Trigger handle

DE: Sprühdüse
 DK: Spraydyse
 EE: Otsik
 ES: Boquilla de pulverización
 FI: Suutin
 FR: Buse de pulvérisation
 GR: Ρύγχος Ψεκασμού
 HU: Spray fej
 IT: Ugello
 NL: Spray opening
 NO: Spraydyse
 PL: Natrysk szczelinowy
 PT: Abertura do difusor
 SE: Spray munstycke
 RO: Aplicator tip spray
 UK: Spray aperture

DE: Überdruckventil
 DK: Trykregulator
 EE: Väljalaskekklapp
 ES: Regulador de presión
 FI: Ylipaineventiili
 FR: Soupape de sécurité
 GR: Βαλβίδα αποσυμπίεσης
 HU: Szelep
 IT: Valvola di sicurezza della pressione
 NL: Overdrukventil
 NO: Trykregulator
 PL: Zawór upustowy
 PT: Válvula de Segurança
 SE: Tryckregulator
 RO: Supapă de presiune
 UK: Pressure relief valve

DE: Überwurfmutter
 DK: Fingerskrue
 EE: Mutter
 ES: Tubo roscado
 FI: Kärkiosa
 FR: Écrou fileté
 GR: Κοντό παξιμάδι
 HU: Menetes anya
 IT: Dado zigrinato
 NL: Kartelmoer
 NO: Fingerskrue
 PL: Nakrętka
 PT: Sede com rosca
 SE: Lasmutter
 RO: Piulița zimțată
 UK: Knurled nut

DE: Behälter
 DK: Flaske
 EE: Balloon
 ES: Recipiente
 FI: Säiliö
 FR: Bouteille
 GR: Φιάλη
 HU: Tartály
 IT: Vontenitore
 NL: Cilinder
 NO: Flaske
 PL: Pojemnik
 PT: Reservatório
 SE: Flaska
 RO: Unitate izolată
 UK: Bottle

$T_{max} = 50 \text{ deg. C}$ $T_{min} = -196 \text{ deg. C}$ $PS = 0,76 \text{ bar}$ $V = 0,35 \text{ l (Mini)}$ $V = 0,50 \text{ l (Maxi)}$

Instruction Manual

DE: Gebrauchsanweisung
 EE: Kasutamisujuhend
 FI: Käyttöohje
 GR: Οδηγίες Χρήσης
 IT: Manuale Utente
 NO: Bruksanvisning
 PT: Manual do usuário
 SE: Instruktion

DK: Betjeningsvejledning
 ES: Manual de utilización
 FR: Notice d'utilisation
 HU: Használati útmutató
 NL: Handleiding
 PL: Instrukcja obsługi
 RO: Manual de utilizare
 UK: Instruction Manual

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UK – CRYOPRO ® INSTRUCTION MANUAL

FILLING THE UNIT

The CryoPro® unit is to be filled with *liquid nitrogen only*. Fill the unit by slowly pouring liquid nitrogen into the top or by using a standard low pressure withdrawal device from a storage tank (Dewar). Depending on the number of lesions to be treated and the duration of the individual treatments the unit should be filled to 25 - 90 % of its maximum capacity. After filling a warm unit, the liquid will boil violently for approx. 30 secs. Hereafter, the top cover can easily be screwed on and tightened firmly. A weak but audible hissing sound from the pressure relief valve indicates that the unit has reached working pressure. This is normal and does not indicate an error condition.

Caution should be taken when refilling a CryoPro® containing residual fluid as the unit is pressurized at 0.74 bar. Slowly unscrewing the top cover will allow the unit to depressurize before the top comes off. This is accompanied by an audible hissing sound until it is fully depressurized.

INTENDED USE

Verucca vulgaris, basal cell carcinoma, cervical cell changes.

OPERATION

The CryoPro® is supplied with a set of standard spray apertures. Optional closed probes are available for contact freezing. All probes and spray apertures are screwed into the permanently affixed knurled nut with fingertight firmness. Do not apply excessive force.

Do never use the CryoPro® without either a spray aperture or closed probe in place. The excessive amount of liquid put out if the unit is used without spray or probe attachment will cause severe damage to the patient.

Spray freezing (REF OS A, B, C, D, SS, BS, SOFT, LL)

The spray apertures include: A (1 mm opening), B (0.75 mm opening), C (0.55 mm opening), D (0.45 mm opening) as well as Bent Spray aperture (BS, 0.55 mm opening) and Straight Spray (SS, 0.55 mm opening) both for areas difficult to reach like behind the ears, lip and outer nose (avoiding gas into the airways).

Select a spray aperture providing an LN₂ output appropriate for the lesion to be treated. In-depth freezing is best obtained when the spray aperture is close to the lesion (5 - 10 mm distance) as the liquid content of the spray is higher close to the aperture. When possible, raised warts should be approached tangentially from at least two sides. This allows the ice ball to move down through the wart while sparing the surrounding tissue. Spraying directly at the lesion (perpendicular to the surface) will cause more lateral spread with less penetration. To avoid this apply the spray intermittently or use a smaller spray aperture.

The SOFT spray aperture is for freezing larger areas at greater distance (up to approx. 5 cm) with higher vapor content. Move the aperture slowly back and forth thus "painting" the surface with vaporized nitrogen.

A Luerlock adapter (LL) facilitates the use of a needle as a spray aperture. The diameter of the needle determines the LN₂ output according to the size of the lesion. Ref. enclosed E280010x.

Contact freezing (REF CP SP, CP 1mm-30mm, CX DS19, CX DS25, CX 19x5, CX 19x25, CX 25x5)

Closed contact probes provide in-depth freezing with minimal lateral spread. Select a probe size fitting the size of the lesion. Make the lesion wet with a drop of water or contact gel prior to freezing to facilitate the thermal transfer between skin and probe tip. Apply the contact probe when still warm applying light pressure and pull the trigger handle. Keep freezing until a 1 mm halo is seen around the probe tip. Then stop the freezing and allow a few seconds for the probe tip to come loose.

Ref. enclosed E280010x, Z111030x.

Freezing times vary according to the aperture or probe used, the size and type of the lesion and the freezing distance (spray only). The user is encouraged to familiarize himself with available literature on cryosurgery using liquid nitrogen.

AFTER USE

At the end of an office day it is highly recommended to store the CryoPro® with the lid on in order to prevent condensation inside the bottle. Build-up of condensation may lead to clogging of the valve assembly.

MAINTENANCE

No preventive maintenance is needed.

DECONTAMINATION

Normally, spray apertures need no cleaning as they are not in contact with the patient. However, should the need for decontamination arise, the same procedure as for closed contact probes applies.

Closed contact probes are in direct skin contact and should be cleaned and steam autoclaved at 121 deg. C for 15 min. after use.

Prevent any liquid or steam from entering the interior of the contact probe via the tube fitting and screw adapter.

Ref. enclosed Z110030x.

STORING LIQUID NITROGEN

To ensure a clean supply of liquid nitrogen the storage tank (dewar) must be completely emptied prior to refilling 3 - 4 times per year. If necessary, water from icocrystals remaining inside the tank can be dried out by creating air circulation (e.g. using a vacuum cleaner for approx. 30 minutes).

WARNINGS

- ◆ Do never use the CryoPro® unit without a spray aperture or closed contact probe in place.
- ◆ After filling, the top cover must be screwed on and tightened securely.
- ◆ The CryoPro® is supposed to be held in a vertical position when used. Liquid nitrogen may escape through the pressure relief valve if the unit is held in a position allowing the liquid inside the bottle to be in contact with the lid.
- ◆ Be careful when opening a pressurized unit. Unscrew the top *slowly*. The unit may still be pressurized even if with no residual liquid left.

The CryoPro® is a surgical instrument. When not in use it must be kept out of reach of children and unauthorized personnel.

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