UK - Decontamination of equipment in direct skin contact

With reference to MHRA MAC manual Part 1 Principles 3rd edition May 2010 the infection risk associated with decontamination of equipment provided by Cortex Technology, and which is intended to be in direct skin contact, is classified as follows:

Risk	Application of item	Recommendation
High	<ul style="list-style-type: none"> in close contact with a break in the skin or mucous membrane 	Cleaning (1) prior to Sterilization (2).

1. Procedure for Cleaning:

Cleaning agent: *Micro 10® Enzyme* by USF Healthcare S.A., Rue François Perréard 4, CH-1225 Chêne Bourg, Geneva, Switzerland, Tél + 41 22 839 79 00, Fax + 41 22 839 79 10, info@usfhc.com.

Follow the manufacturers "Protocol for use" as applicable.

Please observe:

Cryosurgical contact probes and cervical probes shall not be fully immersed into the cleaning solution. Only the part of the probe, which is in direct skin contact during the treatment procedure, shall be fully immersed. Other parts of the item shall be thoroughly cleaned with a cloth soaked in the solution. Any liquid entering the interior of the probe will cause malfunction (blockage) as it freezes up once the probe is used on an LN2 treatment unit.

Other items, where liquid cannot enter the interior parts, may be fully immersed during cleaning.

2. Procedure for Sterilization:

Replace the silicone vent tube on closed contact probes and cervical probes with the protective silicone cap.

Also, attach the silicone stud in the brass thread.

Steam autoclave at 134 deg. C and 2 bar pressure for 3 - 3,5 minutes, 10 minutes drying time.

Remount the silicone vent tube on contact probes and cervical probes prior to use.

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